

Authorized
dealer:



2004-04-20, Z:\Products Public\Loop Amplifiers - Systems\UniVox 2A\English\Series2A_DuG014.pmd

Warning:

Unit may become hot due to excessive input.
Extremely high output level can damage eardrums.
Place the loop pad cord/loop wire cord where you won't stumble over it.
Use in cars may diminish you ability to hear horns & sirens

Bo EDIN AB, specialist in wireless loop systems and assistive devices since 1965.

Bo EDIN AB

Stockby Hantverksby 3
181 75 LIDINGÖ

Office address:

Förrådsvägen 2 B
SE-181 41 LIDINGÖ

Phone:

+46-87 67 18 18

Fax:

+46-87 67 18 20

E-mail:

info@edin.se

Website:

www.edin.se



2004-04-20

UniVox® 2A
Loop Amplifier
Users Guide

Start to check the parts included in the package. How to setup UniVox® 2A is described on page 2 & 3.


UniVox® 2A - Overview

The following parts are included in the UniVox® 2A package:


- A. Amplifier
- B. Power adapter
- C. Loop pad
- D. Input connector cord

UniVox® 2A - Description

Front




Back



UniVox® 2A Placement


Place on side, top, below, or near TV/ Radio. Use velcroadhesive supplied, if desired.



⊗ ⊗ ⚙ UNIVOX 2A + 📺 ⊗

TV / Radio surface

When using a microphone to input sound into Univox, place microphone in front of TV/Radio where sound comes out.



UniVox® 2A Connecting

- A. Connect the power adapter to jack (6)
- B. Connect any earphone to the 3.5-mm stereo jack (7)
- C. Connect the loop pad to jack (8), or
- D. Connect the loop wire to the spring-clip terminals (9). Place the wire around the area you want to cover.
- E. Set the switch (10) for proper setting for LINE (Radio/TV) or MIC (11)
- F. Connect Radio/TV (Line) or microphone to the jack (11)
- G. Connect microphone, if used, to jack (12)

UniVox® 2A Usage / Adjustment

Set your hearing instrument switch to position “T” or “MT”.

Adjust the volume control to desired level. If the sound is distorted the volume is too high. Adjust the bass and treble controls only if necessary to achieve increased clarity.

Volume adjustment (3)

Bass adjustment (1)

Treble adjustment (2)

Loop current indicator (4). A green light indicates the loop/loop pad is working.

Main power indicator (5). A yellow light indicates that the power adapter is properly connected to the amplifier.

UniVox® 2A Trouble shooting

If you experience difficulties in using the system please read these items:

“ON” (yellow) lamp does not light: Check that the adapter is connected both to the wall outlet and to the UniVox.

“LOOP” (green) lamp not glowing: Check that TV/Radio is connected to Univox and loop pad, or wire, is connected.

“LOOP” (green) lamp glowing with no sound: Be sure you are within the wire loop (or over loop pad) and the telecoils in your hearing aids are turned on.

Bad sound quality: Adjust the treble and bass (normally enhancement of bass or treble is not needed).